



## Agar gel immunodiffusion (AGID) test for the detection of EIA anti-p26 antibodies in equine serum and plasma.

- **Clear, bright lines** for easy reading
- **OIE-prescribed protocol**
- May be purchased **with or without agar gel** (the IDvet agar gel allows for clearer precipitation lines)
- **24-month shelf-life**
- **Generates a clean inflexion with OIE reference serum**

### SPECIFICATIONS

<b>Method</b>	AGID (Coggins-test)
<b>Species</b>	Equine
<b>Specimens</b>	Serum and plasma
<b>Antigen</b>	EIAV p26

### ORDERING INFORMATION

<b>Product code</b>	EIA-AGID	EIA-AGID-NOGEL
<b>Kit format</b>	Includes gel	Gel not included
<b>Reactions</b>	200	200

### REFERENCES

1. OIE Terrestrial Manual 2013; 2.5.6; **Equine Infectious Anemia**

### ASSOCIATED PRODUCTS

#### Freeze-dried EIA-AGID reference serum

<b>Product code</b>	MRI-EIA-AGID
<b>Format</b>	0.5 ml vial (freeze-dried)
<b>Description</b>	Freeze-dried EIA positive horse serum to be used as internal

