

Capture ELISA for the detection of the Influenza A nucleoprotein (NP).

- **High sensitivity and specificity**
- **Cost-effective** for the analysis of large sample numbers
- **Applicable to multiple species**, including horses, birds and swine. Please contact IDvet for more information
- **Efficient and robust** on complex samples (faeces and cloacae swabs) which can generate false positive results with immunochromatographic tests
- **Easy-to-use and rapid**, with results in under 1 hour

SPECIFICATIONS

Method	Sandwich ELISA
Species	Multiple species, including birds, pigs, horses and others
Specimens	- Broncho-pulmonary liquids from pigs, horses and birds; - Bird cloacae swab samples; - Bird faeces.
Coated antibody	Anti-nucleoprotein (NP) monoclonal antibody
Conjugate	Anti-NP-HRP conjugate (concentrated 10X)

ORDERING INFORMATION

Product code	INFLAG-2P
Kit format	2 plates
Reactions	192
Plate format	96 individual wells

ASSOCIATED PRODUCTS

Ready-to-use chicken negative serum

Product code	MRINEG-BIRD-RTU
Format	5 ml vial

Description

Ready to use specific pathogen free (SPF) chicken serum pool.

Does not contain any antibodies against any avian pathogen. For use with all the IDvet poultry ELISA kits or for RSA testing.