



# INFECTIOUS BRONCHITIS



## ID Screen® Infectious Bronchitis Competition



Competitive ELISA kit for  
the detection of anti-IBV NP  
antibodies in chicken serum  
and egg yolk samples

**ID.vet**  
Innovative Diagnostics

With you at every step

[www.id-vet.com](http://www.id-vet.com)

# INFECTIONOUS BRONCHITIS

## ID Screen® Infectious Bronchitis Competition

Competitive ELISA kit for the detection of anti-IBV NP antibodies in chicken serum and egg yolk.

- High sensitivity and specificity
- Protocols for serum and egg yolk
- Earlier detection of infected animals
- The stop solution no longer contains H<sub>2</sub>SO<sub>4</sub>, meaning that it is less corrosive

## Specifications

Method	Competitive ELISA
Species	Chickens
Specimens	Serum or egg yolk
Coated antigen	IBV purified antigen
Conjugate	Anti-IBV-HRP conjugate (concentrated 10X)

## Ordering information

Product code	IBVC-2P
Kit format	2 plates
Reactions	192
Plate format	12 x 8-well strips

## Associated Products

### IBV positive freeze-dried serum

Product Code	MRI-IBV
Format	0.25 ml vial (freeze-dried)
Description	Freeze-dried positive chicken serum to be used as internal reference material for quality control.

### SPF negative freeze-dried serum

Product Code	MRINEG-BIRD
Format	0.5 ml vial (freeze-dried)
Description	Freeze-dried SPF negative chicken serum to be used as internal reference material for quality control.

