

## Sandwich ELISA for the detection of p27 antigen of the Avian Leukosis Virus in chicken light albumen or cloacal swabs.

- **Improved analytical sensitivity** with respect to other commercial ELISAS.
- **A wide spectrum of detection**, given that the kit may be used to detect all virus subgroups, including the new K strain circulating in China.

### SPECIFICATIONS

<b>Method</b>	Sandwich ELISA
<b>Species</b>	Chicken
<b>Specimens</b>	Light albumin or cloacal swabs
<b>Coated antibody</b>	Anti-p27 monoclonal antibody
<b>Conjugate</b>	Anti-p27-HRP conjugate (concentrated 10X)

### ORDERING INFORMATION

<b>Product code</b>	ALVAG-5P	ALVAG-10P
<b>Kit format</b>	5 plates	10 plates
<b>Reactions</b>	480	960
<b>Plate format</b>	1 x 96-well whole plates	1 x 96-well whole plates